

Order at your participating dealer today:
www.renfert.com/dealers

Renfert

Basic master

Fine sandblasting unit with 2 tanks

Precise blasting with the special mixing chamber technology (Venturi principle) ■ Large volume blasting chamber provides maximum freedom of movement ■ Ventilation slots ensure optimal ventilation of the blasting chamber and a clear view (PerfectView-Concept)



Basic master,
1 x 25 – 70 µm,
1 x 70 – 250 µm,
Art. No. 2948 2025

€ 1.040,–

Basic quattro IS

Sandblasting unit
with quick stop technology

Controlled and exact sandblasting of finely detailed objects thanks to integrated quick stop technology 'IS' (Immediate Stop) ■ Precise and reliable processing at best possible cost effectiveness with the special mixing chamber technology



Basic quattro IS,
1 x 25 – 70 µm,
1 x 70 – 250 µm,
Art. No. 2959 0000

€ 1.640,–

Vario jet + Vortex compact 3L

Recyclable sandblaster
+ Wet and dry extraction unit

Vario jet: Exact coordination of the essential sandblasting parameters (pressure, time, grit size) ■ Intelligent extraction system: only unusable fine dust is extracted ■ High degree of freedom of movement due to easily removable sandblasting basket ■ **Vortex compact 3L:** Optimum protection against fine dusts (99.9% filtration efficiency) using a 2-stage filter system with M Class fine filters, optional also H/Hepa filters

Vario jet,
Art. No. 2961 0000
+ Vortex compact 3L,
Art. No. 2924 5000

€ 3.300,–



ΟΙ ΤΙΜΕΣ ΔΕΝ ΠΕΡΙΛΑΜΒΑΝΟΥΝ ΦΠΑ