

A-dec 300 Features/Benefits

Full Radius capability

Pedestal Design

Lower retail price

Improved aesthetic

A-dec 300 ancillary integration

Back actuator disconnect

Integrated USB cable

Manual height adjust arm

Easy amalgam separator access

2L QD water bottle

Hybrid drive

Floating equipment

Dual Radius front mount

Improved comfort

Trendelenburg

Improved stability

Optional assistants arm touchpad

Optional two position armrests

Fully upholstered toeboard

Quick support center access

Platform design

new A-dec 311 chair

- fully upholstered toe board for greater comfort and aesthetic appeal
- optional armrest with simple forward rotation
- virtual pivot for improved comfort



A-dec Module Compatibility

311 Radius System

Back B

Products that are connected to the back-mount of the chair.



461 support center/cuspidor



551 assistant's instrumentation (shown with 461, compatible with or without cuspidor)



351 assistant's instrumentation (not compatible with cuspidor)



Front F

Products that are connected to the front-mount of the chair.

Dual Mount Choices



332 Traditional delivery system (shown)



333 Continental delivery system



336 Traditional delivery system



572L dental light (shown)



572 dental light



382 monitor mount

Light post monitor options



381 monitor mount (mounts to light post)



482 monitor mount (shown, mounts to light post)

Dedicated Left-Right Post-Mount System (for any dental chair accepting a post-mount system)

A-dec Module Compatibility

Side S

Products that are connected to the side-mount of the chair. Left or right option is available.



362 post-mounted support center with cuspidor or without cuspidor (shown)



334 Traditional delivery system (shown)



335 Continental delivery system



352 telescoping assistant's instrumentation (shown)



353 cuspidor-mounted assistant's instrumentation



571L dental light



571 dental light (shown)



381 monitor mount (mounts to light post)



482 monitor mount (mounts to light post)



(Shown on an A-dec 311 chair)

A-dec Module Compatibility

311 Dedicated Left-Right Pedestal System – International only

Side S

Products that are connected to the side-mount of the chair. Left or right option is available.



363 post-mounted support center with cuspidor (shown) or without cuspidor



334 Traditional delivery system



335 Continental delivery system (shown)



352 telescoping assistant's instrumentation



353 cuspidor-mounted assistant's instrumentation (shown)



571L dental light



571 dental light (shown)



381 monitor mount (mounts to light post) (shown)




482 monitor mount (mounts to light post)




Choice & Compatibility - 300 Dedicated

► Dedicated Left- or Right-Side Mounted Modules



Choice & Compatibility - 300 Radius



A-dec Module Compatibility

511 Radius System



Back B

Products that are connected to the back-mount of the chair.



461 support center/cuspidor



561 support center/cuspidor
(compatible with or without cuspidor)



551 assistant's instrumentation (shown)



351 assistant's instrumentation
(not compatible with cuspidor)



571L dental light



571 dental light (shown)



581 monitor mount (shown)

Front F

Products that are connected to the front-mount of the chair.

Single Mount Choices



532 Traditional delivery system



533 Continental delivery system
(shown)

Dual Mount Choices



332 Traditional delivery system



333 Continental delivery system



336 Traditional delivery system



382 monitor mount



572L dental light



572 dental light

Light Post Monitor Options (Back or Front)



381 monitor mount



482 monitor mount (shown)

► Upholstery – sewn and seamless.

Sewn Upholstery



Seamless Upholstery

