

**E U R U S**  
**Options & Accessories**

Chair



Double articulating headrest



Power headrest

Dentist Element



Holder type



Rod type



Place type



Belmont handpiece



Syringe Belmont BT14 & 77



Wireless Foot controller

Unit Body



Vacuum handpiece & Saliva ejector Belmont & Durr



Twin Vacuum(Two solid correctors)



Sensored Cupfiller



Clean water system with 2 liter bottle



LED indicator



Separated Junction Box



Monitor bracket (not incl.monitor)



Vacuum line clean system (Pedestal type only)





Flushing system (Integrated type)

Light



Patient mirror



Stool



Stool for Dentist

● Standard | ○ Selectable | - Not available | *PD* Pedestal type only

| Chair                                                                                 | S1 & S3 Contour | S4 Folding legrest |
|---------------------------------------------------------------------------------------|-----------------|--------------------|
| Double articulating headrest                                                          | ○               | ○                  |
| Power Headrest                                                                        | ○               | ○                  |
| Armrest at left side only (Fixed)                                                     | ○               | -                  |
| Armrest at both side, Left side fixed & Right side rotated                            | ○               | ●                  |
| Integrated stick switch (2 Controls, Manual & Auto mode chair switch)                 | ○               | -                  |
| Integrated stick switch (3 Controls, Manual, Auto mode chair switch & Power headrest) | ○               | ●                  |
| Upholstery Seam type                                                                  | ○               | ○                  |
| Upholstery Luxury Seam type                                                           | ○               | ○                  |
| Upholstery Seamless type                                                              | ○               | ○                  |

| Dentist Element                                                                  | Holder/Place | Rod |
|----------------------------------------------------------------------------------|--------------|-----|
| Touch panel screen, 5 inch                                                       | ●            | ●   |
| Large table with silicon mat                                                     | ●            | -   |
| Small size tray holder (size: 290x190mm)                                         | -            | ○   |
| Large size tray holder (size: 380x290mm)                                         | -            | ○   |
| 3-way syringe                                                                    | ○            | ○   |
| 6-way syringe                                                                    | ○            | ○   |
| Midwest 4-way turbine hose                                                       | ○            | ○   |
| Midwest 4-way turbine hose with fiber optic                                      | ○            | ○   |
| Micro motor (Air/Electric)                                                       | ○            | ○   |
| Scaler (Air/Electric)                                                            | ○            | ○   |
| Water heater for handpiece                                                       | ○            | ○   |
| Foot controller (wired)                                                          | ○            | ○   |
| Foot controller (wireless)                                                       | ○            | ○   |
| Intra oral camera                                                                | ○            | ○   |
| USB Port                                                                         | ○            | ○   |
| Stainless Waste Receptacle/ Cotton Container with Stand /Vials x 5pcs with Stand | ○            | ○   |
| Intra oral Film viewer                                                           | ○            | ○   |
| Panorama Film viewer                                                             | ○            | ○   |
| Curing light                                                                     | ○            | ○   |

| Assistant Element                 | Holder/Place | Rod |
|-----------------------------------|--------------|-----|
| Standard type Assistant arm       | ○            | ○   |
| Ratchet gear type Assistant arm   | ○            | ○   |
| 3-way syringe                     | ○            | ○   |
| Vacuum & Saliva ejector handpiece | ○            | ○   |
| Twin Vacuum                       | ○            | ○   |

| Unit Body                                   | Holder/Place  | Rod           |
|---------------------------------------------|---------------|---------------|
| Electric timed Cupfiller & Bowl flush       | ○             | ○             |
| Sensored Cupfiller & Bowl flush             | ○             | ○             |
| LED indicator                               | ○             | ○             |
| Rotating cuspidor bowl, porcelain           | ●             | ●             |
| Water heater 24V                            | ○             | ○             |
| Water disinfection system                   | ○             | ○             |
| Clean water system with 2 liters water tank | ○             | ○             |
| Separator                                   | ○             | ○             |
| Amalgam separator                           | ○             | ○             |
| Vacuum line clean system                    | ○ / <i>PD</i> | ○ / <i>PD</i> |
| Flushing system without flushing tray       | ●             | ●             |
| Flushing system with flushing tray          | ○             | ○             |
| Flushing system (Integrated type)           | ○             | ○             |
| Monitor bracket                             | ○             | ○             |
| Service outlet                              | ●             | ●             |
| Connector for Water                         | ○             | ○             |
| Service outlet/Connector for Air            | ○             | ○             |
| Separated Junction Box                      | ○             | ○             |

| Light          |   |
|----------------|---|
| Eurus Light    | ● |
| Patient mirror | ○ |

| Stool                     |   |
|---------------------------|---|
| Eurus Stool for Dentist   | ○ |
| Eurus Stool for Assistant | ○ |

## Dimensions

**S1**



**S3**



**S4**

