

OPTOMIC®

producing quality

Introducing the new
extension system for
OP-DENT

**Optimal ergonomics
for your microscope**




A 45-DEGREE EXTENSION
INCREASES THE DISTANCE
BETWEEN YOUR
EYES AND THE
PATIENT'S MOUTH.


THIS ALLOWS YOUR
SHOULDERS AND ARMS
TO WORK AT CORRECT
ANGULATION



producing quality



25° righth rotation



25° left rotation

THE EXTENSION SYSTEM PROVIDES **50-DEGREE ROTATION**.
YOU CAN **MAINTAIN YOUR EYES IN THE HORIZONTAL**
WHILE **FOCUSING** THE MICROSCOPE
ON THE **CORRECT WORK AREA**.



producing quality



producing quality

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain

Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14

export@optomic.com • www.optomic.com

